

Microbiology [Glassware & Plasticware]

Sr]	Name of the Test	Make	Requirement
01	Flask [1000ml glass autoclavable	Borosil / Borosilicate	5 Nos
02	Flask [500ml glass autoclavable	Borosil / Borosilicate	5 Nos
03	Flask [250ml glass autoclavable	Borosil / Borosilicate	5 Nos
04	Flask [100ml glass autoclavable	Borosil / Borosilicate	5 Nos
05	Beaker [500 ml] glass autoclavable	Borosil / Borosilicate	2 Nos
06	Beaker [150 ml] glass autoclavable	Borosil / Borosilicate	5 Nos
07	Beaker [100 ml] glass autoclavable	Borosil / Borosilicate	2 Nos
08	Screw Cap bottle [1000ml] glass autoclavable	Borosil / Borosilicate	5 Nos
09	Measuring Cylinder [100ml] glass autoclavable	Borosil / Borosilicate	2 Nos
10	Test tube [18 x 150 ml] glass autoclavable	Borosil / Borosilicate	2500 Nos
11	Test tube [15 x 125 ml] glass autoclavable	Borosil / Borosilicate	2500 Nos
12	Test tube [12 x 100 ml] glass autoclavable	Borosil / Borosilicate	2500 Nos
13	Blood culture bottle 60 ml with rubber seal & graduated Glass autoclavable	Borosil / Borosilicate	2500 each
14	Blood culture bottle 60 ml with rubber seal & graduated Glass autoclavable	Borosil / Borosilicate	6000 Nos
15	Pipette [2 ml] glass autoclavable	Borosil / Borosilicate	5 Box
16	Petri dish [100 mm]	Laxbro	1 Nos
17	Vial Screw Cap 17mm x 75mm [1 x 500]	Hi-media	12000 Nos
18	Autoclavable Plastic Spetum Container Universal [30ml]	----	20000 Nos
19	Sterile Culture swab wooden stick [Hi-media]	PW 005	35000 Nos
20	Cover slip [22 x 22mm] glass	Blue Star	400 Pkt
21	Staining bottle glass autoclavable	----	15 Nos
22	Micro Pipette tips 20-200 ul	Tarson	10000 Nos
23	Micro Pipette tips 1-10 ul	Tarson	1000 Nos
24	PCR Tube with flat lid 0.2 ml capacity	Tarson	4000 Nos
25	Spectrophotometer	-----	250 Nos
26	Test tube rack [Steel/Aluminium] 125 x 15 mm	-----	4 Nos
27	Test tube rack [Steel / Aluminium] 100 x 12 mm	-----	4 Nos
28	Nicrome loop 4 [4mm diameter,double wound Caliberated to 0.01ml Hi-media LA 019	Hi-media LA040	100 Nos [10 x10]
29	Metal loop holder [w/o loop] brass rod holder with heat Resistant handle where any size of nichrome wire can be Inserted	Hi-media LA040 [1 x 24 Nos]	48 loop holder [2 x 24]
30	Parafilm Parafilm taps to seal petridish [para film]	LA 017 INO	2 x 1 No